# BEAUTIFLY

White

# Smile

series

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Power:	1.8W
Voltage:	5V
Waterproof:	IPX7
Vibration:	38000/min
Battery capacity:	1200mA

### Safety

Read the instructions carefully before use. Follow these safety guidelines when using the device.

### Description of the device [page 2 - f1]

1. bristle
2. brush head
3. power button
4. charge indicator
5. gear display
6. USB charging cable

### Using method [page 2 - f2]

1. Insert the brush head into the output shaft, leaving a slight gap between the brush head and the handle. This allows the brush head to vibrate properly.
2. Take a small amount of toothpaste to prevent it from splashing.
3. Wet the toothbrush head before putting it in the mouth, then press the power button.

Choose the gears according to your needs. The toothbrush automatically switches to cleaning mode.

4. Place the toothbrush bristles at a 45° angle to your teeth, so they are in line with the gum line. With circular motions, gently move the brush head. Do not press or scrub hard.

5. You will experience a brief pause every 30 seconds during use. It reminds you to move to another part of your mouth to achieve a more effective cleaning effect. After two minutes of brushing, the toothbrush automatically switches off.

6. Rinse your mouth with clear water after cleaning.

### Brush operating modes [page 2 - f3]

1. Cleaning
2. Polishing
3. Gentle
4. Whitening
5. Gum Care

### Battery level indicator light

Red light flashes when it requires charging.

### Over-discharge protection

The output will be stopped when the voltage reaches 2.8V.

### Timing function

Regardless of the mode of operation, the working time is two minutes.

### Cleaning and maintenance [page 2 - f4]

#### 1. Brush Head

After each use, clean the brush head and bristles and let them air dry naturally. Do not pull the bristles; otherwise they will deform and drop off.

#### 2. Main Body

The main body of the toothbrush can be washed with water and dried naturally. Do not immerse it in water. Do not use water that is more than 40°C or an acid cleaning solution.

#### 3. Charging Base

Unplug the charger and wipe the entire surface of the charging base with a damp cloth.

#### 4. Brush head replacement

Replace the brush head every three months or sooner if it becomes worn.

### Note

1. Before first use, start the brush up. If electricity remains, use it up before charging.
2. Charge the toothbrush for more than 12 hours if you do not use it for 6 months. After fully charged, it can be used for about 30 days.

## Warnings

- Do not pull the charger above or near the water container, such as the bathtub, pool, etc.
- Do not use the device if the electric cord or power supply has been damaged or is working improperly.
- If the toothbrush is not used for a long period of time, the indicator light will not appear for a short time during charge. Wait for 6 or 10 minutes and the light will go back to normal.
- When charging, the contact between the teeth and the body of the device is warmer, which is normal. If the heat emitted from this area is high and you feel discomfort, stop using the device.
- This product requires a matched brush head and a matched charger. Cut off the power supply when the charging stand is not in use or before cleaning.
- This product is intended for cleaning teeth only. Do not use it for other purposes.
- This product is only for personal use. It can't be used by dental clinics or institutions for more than one patient.
- Do not expose the toothbrush to direct sunlight. Place it in a dry and ventilated environment.
- Children and people with limited physical, sensory, or mental abilities are allowed to use this appliance, provided they have been given supervision or instruction from a person responsible for their safety.
- Children should be supervised to ensure they do not play with the appliance.
- Don't let hair care products (such as gelatin, oil) or hand sanitizers contact with this product so as not to damage the surface.
- Clean the power plug regularly to keep dust and dirt from building up and reducing the effect of lemon elimination, which can lead to overheating or electric shock and unnecessary losses.
- When scrapping the main body of the toothbrush, follow the local regulations on scrapping and recycling of small household appliances or take back the manufacturer's unified treatment.

## Storage

- Keep it dry and clean.
- Keep it away from high temperatures.
- Keep it out of reach of children.
- If you do not want to use this product for a long time, you should unplug the power plug.
- Place it in a dry, ventilated place. Do not expose the toothbrush to direct sunlight.

## More information about the company and products [page 2 - f5]

- [www.beautify.eu](http://www.beautify.eu)
- [www.youtube.com/beautify.eu](http://www.youtube.com/beautify.eu)

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**[EN] Rules for handling electrical and electronic equipment:** Electrical and electronic equipment as well as batteries and accumulators must not be combined with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Dispose of waste electrical and electronic equipment as well as batteries and accumulators at your local separate collection point or to your dealer. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Incorrect handling or damage to them may lead to damage to health or environmental pollution when the equipment is further disposed of or recycled. The device was placed on the market after August 13, 2005

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**[EN] Summary declaration of conformity:** The manufacturer hereby declares under its sole responsibility that the above product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances.

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