# BEAUTIFLY

Derma  
Blush  
series

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[beautifly.eu](https://www.beautifly.eu)

[youtube.com/beautiflyeu](https://www.youtube.com/beautiflyeu)

Tension:	5V
Current:	450mAh
Battery:	Li-ion battery
Power:	1.6W
Size:	163.5x52x38mm
Libra:	205g
Suction power:	50-60 kPa

### **Part description - device (page 2, figure f1)**

1. Charging indicator
2. Suction power indicator
3. Power button/suction power adjustment
4. USB port

### **Parts description - heads (page 2, figure f2)**

1. **Small round suction head.** Low suction power For people who have dry or sensitive skin and those who start treatments with the device.
2. **Oval suction head.** Medium suction power/special shape. For lifting and firming treatments.
3. **Round suction head.** Gentle/medium suction power. For people who have oily skin, acne or blackheads.
4. **Small microcrystalline suction head**
5. **Oval microcrystalline suction head**
6. **Round microcrystalline suction head**

### **Service**

#### **On/Off**

Press and hold the power button for 3 seconds, the product is turned on, the display is on, and the suction power displays 1 level (soft).

### **Power regulation**

Short press the power button to adjust the suction power. Each short press increases the suction power by one level, in the order: soft > normal > strong > ultra strong.

### **Application**

- Discoloration
- Scars
- Stretch marks
- Enlarged, highly visible pores
- Blackheads and acne (avoid inflammatory and purulent lesions)
- Wrinkles
- Skin lacking a healthy color with the first signs of photoaging

## Contraindications

- Active purulent acne
- Rosacea
- Retinoid therapy
- Acid treatments, lasers - at least 4 weeks before microdermabrasion
- Psoriasis
- Viral, bacterial and fungal diseases of the skin
- Vascular skin
- Tendency to scarring
- Broken skin (fresh wounds)
- Skin tumors or suspicious moles in the treated area
- Fillers in the face - at least 4 weeks before performing microdermabrasion, otherwise the preparation may move.
- Do not wash the device with water, the device is not waterproof.

Microdermabrasion (page 2 - figure f3) - move in the direction of the arrows

Lifting massage (page 2 - figure f4) - move in the direction of the arrows

Set the suction head to a 90 degree position and slowly move it over your face. Do not hold in the same part for a long time. (Before removing the suction head from the skin, tilt it to the side to release the pressurized air. This will allow you to safely remove it without irritating the skin. Do not pull or move it directly to avoid pulling the skin and causing redness.)

A detailed description can be found at:

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EN: Electrical and electronic equipment, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out wheeled bin symbol on the packaging. Used electrical and electronic equipment, as well as batteries and accumulators, should be taken to local waste collection points or returned to the seller, and further details can be obtained from your municipality. Electrical and electronic equipment, as well as batteries and accumulators, may contain hazardous substances, mixtures and harmful components that are dangerous for the environment and human health. Improper handling or damage to these items can lead to health problems or environmental pollution during the disposal or recycling of equipment.

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**EN:** The manufacturer hereby declares, on their sole responsibility, that the product Derma Blush complies with the Directive of the European Parliament and of the Council 2014/35/EU of 26 February 2014 on the harmonization of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, the Directive of the European Parliament and of the Council 2014/30/EU of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, and the Directive of the European Parliament and of the Council 2011/65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU with regard to the list of restricted substances. The full text of the Declaration of Conformity is available at:

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