# BEAUTIFLY

PowerFit

# Fitness

series

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Battery capacity:	350mAh
Input voltage:	5V
Charging power:	1,5W
Charger:	USB
Full charge time:	1,5h
Pulse frequency:	4-60Hz
Number of modes:	7
Number of intensity levels:	15
Waterproof level:	IPX5
Net/gross weight:	148/390g
Included:	device, USB cable, manual, warranty card

Description of the device parts [see page 2, f1]

1. Button "+" to increase the intensity
2. Button "-" to decrease the intensity
3. Electrode
4. Connection cable
5. Numerical value indicator
6. Charge indicator
7. USB connector
8. ON/OFF and mode button
9. Electrode

Safety precautions

- Please review the following precautions and contraindications before using.

Attention

- The device is used by holding it in your hands, it is not recommended to use the device in any other way or for other purposes.
- Remove jewelry and metal objects from your body before use;
- Stop using immediately if you experience discomfort, pain or any abnormality;
- It is strictly forbidden to wet the device or submerge it in water;
- Do not use in rooms with high humidity (such as bathrooms, saunas, swimming pools, etc.);
- Do not use the device near flammable materials and heat sources;
- Do not leave the device plugged in when not in use. Unplug the device from the power supply immediately after use;
- Keep the device out of the reach of children.
- The device is not intended for use by people with mental, physical or sensory disabilities;
- In the event of a malfunction, do not disassemble, repair or modify the device yourself - failure to do so will void the supplier's warranty to the purchaser;
- The device is not intended for medical purposes and is intended for preventive purposes.

Contraindications

- heart disease;
- the presence of a pacemaker;
- high fever;
- acute disease;
- high or low blood pressure;
- pregnancy, lactation;
- it is forbidden to use in the area of painful skin areas, wounds, abrasions and cuts;
- skin inflammation;
- oncological disease;
- varicose veins;
- tuberculosis;
- wet body.

## Charging

1. Make sure to fully charge the device before the first use. Full charging of the device usually takes 1.5 hours (when it is completely discharged).
  2. To start the charging, connect the device to a power source using a USB cable.
  3. While charging, one charge indicator symbol will blink. When the device is fully charged, the three charge indicator symbols will light up.
  4. The fully charged device can be continuously used for 5 hours.
- \*The device is not usable when being charged.

## Device control

The Light indicators on the display show the selected mode and intensity level.

- EMS mode from 1 to 7
  - EMS intensity level from 1 to 15
1. Press and hold the ON/OFF button for 2 seconds. The device will start on mode 1 and intensity level 1.
  2. Press the ON/OFF button repeatedly to select the desired mode.
  3. To adjust the required intensity of the selected mode, use buttons +/- to the required level.
  4. If the mode is switched while the device is running, the 1st intensity level will be automatically set.
  5. If you release one of the handles while the device is running, after about 3 seconds it will be set to the 1st intensity Level, and then, after about 30 seconds, it will beep twice and the device will automatically shut off.
  6. To turn on the device, press and hold the ON/OFF button for 2 seconds.

## Storage

- Keep the device dry and clean.
- Keep away from high temperatures.
- Keep away from children.
- Make sure the massager is turned off and unplugged.
- Store the device in a cool, dry place.
- Direct sunlight is not recommended.

## Cleaning

- Use a soft cloth slightly moistened with water to clean the surface of the device.
- Do not immerse the device in water.
- Unplug the massager when cleaning.
- Do not use oil, solvents or harsh detergents to clean the massager.

More information about the product is available on the manufacturer's website at: [see page 2, f2]

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EN: Electrical and electronic devices, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out trash can symbol on the packaging. Used electrical and electronic devices, batteries, and accumulators must be taken to local selective waste collection points or to the retailer, and more information can be found in your municipality. Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components harmful to the environment and human health. Improper handling or damage to them can lead to health hazards or environmental

pollution during further disposal or recycling.

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EN: The manufacturer hereby declares, at its sole responsibility, that the PowerFit Fitness product complies with the Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as with the Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU with regard to the list of restricted substances. The full text of the declaration of conformity is available at the following website:

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