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GE GARETT

Beauty Powerful Mini

O produkcie

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Beauty Powerful Mini

About product

Garett Beauty Powerful Mini is a wireless modern vibration massage device that relaxes sore, stiff muscles, improves blood circulation and the general condition of the soft tissues of the body.

Product description

4 replaceable types of massaging heads:

'U' head
for vertebral massage

flat head
for joint/tendon massage

ball head
for muscle massage

precision head
for all parts of the body

Notes for users

1. Before using this product, please read this manual carefully and keep it for future reference.
2. The product is charged with the charging kit provided.
3. It is forbidden to use this product while charging.
4. The time of continuous use of the device should not exceed 1 hour.
5. The product is not waterproof, do not submerge it in water.
6. Do not attempt to disconnect the device from the outlet or use it with wet hands.
7. To avoid accidents, do not touch the metal part of the plug.
8. Do not damage, modify, bend or twist the power supply cord. This may damage it.
9. Avoid getting the charging port and power supply cord dirty.
10. You should stop the charging if the power plug or cord is damaged or it becomes hot during charging.
11. Clean the device to avoid dust and dirt build-up.
12. Do not store the product in the bathroom, humid places, at high temperature, as this may cause electric shock or fire.
13. Do not destroy the product, throw it or modify it. In the event of damage, please stop using it.
14. Do not use the device around irritated, allergic or inflamed skin. In such case, please contact your doctor.
15. Children should not be allowed to use the device.
16. Pregnant women or women during the menstrual cycle should not use the product without consulting a doctor first.
17. Do not use the device under the influence of alcohol or if you feel unwell.

Instructions for use

1. Make sure the device is charged before using it for the first time.
2. If there is no response to the power button, connect the device for charging.
3. When the device is discharged, it takes about 3 hours to fully charge the battery.
4. The first two numbers on the display indicate the battery charge level and inform about active charging.
5. To turn on the device press . The massager will start the operation in the default mode - level 00. Press the power button again to turn off the device.

6. Using the power control buttons, you can adjust the intensity of the massage according to your individual needs.
7. To change the head, grasp the replaceable end and while gently twisting it, pull it out of the opening.
8. When changing the massage head, in order to improve its operation and prevent seizures, you can use a special agent with the consistency of oil, gently rubbing the walls of the hole from which it was taken.
9. Installation is similar: position the head correctly, and then gently push it until it is motionless.
10. To turn off the device, press the button . The switch-off will be signalled by an acoustic signal.

Note: The default device runtime is 15 minutes. If you want to extend the use, you should restart it.

Massage techniques

Massage of the lower extremities muscles:

Anterior muscles of thigh:

- *Type of head: flat
- *Massage angle: 90°
- *Massage direction: from top to bottom
- *Time of massage: 30 - 45 s
- *Number of repetitions: 2

Posterior muscles of thigh:

- *Type of head: flat, precision
- *Massage angle: 90°
- *Massage direction: from top to bottom
- *Time of massage: 20 - 30 s
- *Number of repetitions: 3

Medial muscles of thigh:

- *Type of head: precision
- *Massage angle: 45°
- *Massage direction: from top to bottom
- *Time of massage: 15 - 25 s
- *Number of repetitions: 2

Anterior muscles of tibia:

- *Type of head: precision
- *Massage angle: 90°
- *Massage direction: from top to bottom
- *Time of massage: 25 - 35 s
- *Number of repetitions: 2

Posterior muscles of tibia:

- *Type of head: ball
- *Massage angle: 90°
- *Massage direction: from top to bottom
- *Time of massage: 15 - 30 s
- *Number of repetitions: 3

Anterior/posterior muscles of foot:

- *Type of head: flat, ball, 'U', choice depending on the needs
- *Massage angle: 90°
- *Massage direction: from top to bottom
- *Time of massage: 15 - 30 s
- *Number of repetitions: 3

Massage of the buttock muscles:**Gluteus minimus muscles:**

- *Type of head: flat, ball
- *Massage angle: 90°
- *Massage direction: from inside to outside, from top to bottom
- *Time of massage: 30 s - 1 min
- *Number of repetitions: 4

Gluteus medius muscles:

- *Type of head: flat, ball
- *Massage angle: 90°
- *Massage direction: from inside to outside, from top to bottom
- *Time of massage: 30 s - 1 min
- *Number of repetitions: 4

Gluteus maximus muscles:

- *Type of head: flat, ball
- *Massage angle: 90°
- *Massage direction: from inside to outside, from top to bottom
- *Time of massage: 30 s - 1 min
- *Number of repetitions: 4

Massage of the sacroiliac joint:**Gluteus minimus muscles:**

- *Type of head: flat, ball
- *Massage angle: 90°
- *Massage direction: from inside to outside, from top to bottom
- *Time of massage: 1 - 1,5 min
- *Number of repetitions: 1

The selected techniques described above can be used as exemplary treatment plans. The choice of massage parameters depends on the individual needs of the user.

Health tips

The massage therapy and the correct use of the device can quickly reduce pain caused by joint/muscle strain, neuralgia, muscle spasms and other musculoskeletal disorders. Most soft tissue injuries are related to these injuries or to incorrect posture. Most problems can be treated with the Beauty Powerful Mini massager. **Avoid excessive stress during the treatment, do not use the massager around sensitive areas such as the head, face and cervical vertebrae.** Stop the massage if inflammation, swelling, or increased pain occurs. The vibration massage is a technology that allows you to heal ailments, muscular problems, rest after a hard workout or relax in your free time.

Indications	Contraindications
Pain and spasms caused by muscle injuries, sprains, strains - supporting the treatment and regeneration of tissues.	Swelling, bloodshot or bleeding parts of the body.
Slightly swollen skin areas - supporting the regeneration of tissues.	Heart diseases, use of heart pacemaker.
Sore muscles - decreasing and reducing the build-up of lactic acid.	Pregnant women, children.
Muscle tiredness - increasing the mobility, reduction of pain.	Serious injuries - it is recommended to consult a doctor first.
A feeling of tiredness, pain - relaxing massage.	

Repair and maintenance

1. In the event of damage or failure, please stop using this product and contact the Garrett service.
2. Any attempts to do the repair on your own may damage the device and invalidate the warranty.
3. The device is not waterproof. To avoid its damage please do not wash it with water.
4. If the device is not used for a long time, it is recommended to charge it at least 1-2 times a month and store it in a dry and cool place to extend its life.

Technical specification

Product name: Garrett Beauty Powerful Mini

Battery capacity: 2000 mAh	Vibration freq.: 2000 - 3000 rpm
Charging time: about 3 hours	Rated power: 25 W
Input voltage: 7,4 V/1 A	Rated voltage: 3,7 V
Weight: 395 g	Dimensions: 17,5 x 14,7 x 5 cm

More information on this device is available on the website www.garett.eu. Visit the website www.garett.eu for learning about more products and accessories.

Warranty/Complaint

The product is covered by 24-month manufacturer's warranty. The warranty does not cover: defects and damage caused by the user's fault or as a result of improper use (not in accordance with the intended use), mechanical damage and resulting from the use of accessories other than those included in the set. When complaining about the product, please complete the warranty card attached to the device and send it together with the proof of purchase to the service address provided in the warranty card.

Contact details available on the website:

www.garett.eu

Producer:

GARETT sp. z o.o.
ul. Targowa 18/1413
25-520 Kielce
www.garett.eu

Correct Disposal of this Product
(Waste Electrical & Electronic Equipment)

The marking on the product or in related texts indicates that it should not be disposed of with other household wastes at the end of its working life. To avoid harmful effects on the environment and human health due to uncontrolled waste disposal, please separate the product from other types of waste and recycle responsibly to promote the reuse of material resources as permanent practice. Household users may obtain information on where and how to recycle this product for the environment in a safe way from the retailer where they purchased the product or from local authorities. Business users should contact their supplier and check the terms of the purchase contract. The product should not be disposed of with other commercial waste.

RoHS

Produced in P.R.C.